



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Calico Triangular Bandage Specifications

Catalogue No. /Reference No. /Products Code	JC01010
Material	100% unbleached cotton, muslin cotton or calico
Cotton Yarn	40's, 32's, 21's
Gauze Mesh	40's, 60x58, 50x48, 43x27, 45x25, etc
Size	90cm x 90cm x 127cm (36" x 36" x 51"); 96cm x 96cm x 136cm; 100cm x 100cm x 141cm (40" x 40" x 56"); 106cm x 106cm x 150cm;
Ply	6ply
Sterile or Not	No Sterile
Rules & Classification	Rule 1, Class I





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Calico Triangular Bandage Specifications

Reference No/ Product Code	Description & Size	Size	Unit Package	Outside/Carton Package
JC010101	Calico Triangular Bandage, with different fabric mesh based on the customer request, with two safety pins, 1pc/zipper bag with printing label, 250 pcs/carton; Vacuum Package is available.	90cm x 90cm x 127cm (36" x 36" x 51")	1 Piece/Zipper Bag	250 Pcs/Carton
JC010102		96cm x 96cm x 136cm	1 Piece/Zipper Bag	250 Pcs/Carton
JC010103		100cm x 100cm x 141cm (40" x 40" x 56")	1 Piece/Zipper Bag	250 Pcs/Carton
JC010104		106cm x 106cm x 150cm	1 Piece/Zipper Bag	250 Pcs/Carton